

Reverse and multiple stable isotope probing to study bacterial metabolism and interactions at the single cell level

Yun Wang¹, Yizhi Song², Yifan Tao³, Howbeer Muhamadali⁴, Royston Goodacre⁴, Ning-Yi Zhou⁵, Gail M. Preston⁶, Jian Xu^{1}, Wei E. Huang^{2*}*

1. CAS Key Laboratory of Biofuels and Shandong Key Laboratory of Energy Genetics and Single Cell Center, Qingdao Institute of BioEnergy and Bioprocess Technology, Chinese Academy of Sciences, Qingdao, 266101, P. R. China.
2. Department of Engineering Science, University of Oxford, Parks Road, Oxford, OX1 3PJ, United Kingdom.
3. Department of Operative Dentistry and Endodontics, Guanghua School and Hospital of Stomatology, Sun Yat-sen University, Guangzhou, 510055, China.
4. School of Chemistry, Manchester Institute of Biotechnology, University of Manchester, 131 Princess Street, Manchester, M1 7DN, United Kingdom.
5. State Key Laboratory of Microbial Metabolism and School of Life Sciences & Biotechnology, Shanghai Jiao Tong University, Shanghai, P. R. China
6. Department of Plant Sciences, University of Oxford, South Parks Road, Oxford OX1 3RB, United Kingdom.

*Corresponding author: Jian Xu and Wei E. Huang

xujian@qibebt.ac.cn and wei.huang@eng.ox.ac.uk

Telephone: +44 (0)1865 283786, Fax: +44 (0)1865 374992

ABSTRACT

The interactions between microorganisms driven by substrate metabolism and energy flow are important to shape diversity, abundance, and structure of a microbial community. Single cell technologies are useful tools for dissecting the functions of individual members and their interactions in microbial communities. Here, we developed a novel Raman stable isotope probing (Raman-SIP), which uses Raman micro-spectroscopy coupled with reverse and D₂O co-labelling to study metabolic interactions in a two-species community consisting of *Acinetobacter baylyi* ADP1 and *Escherichia coli*. This Raman-SIP approach is able to detect carbon assimilation and general metabolic activity simultaneously. Taking advantage of Raman shift of single cell Raman spectra (SCRS) mediated by incorporation of stable-isotopic substrates, Raman-SIP with reverse labelling has been applied to detect initially ¹³C-labelled bands of ADP1 SCRS reverting back to ¹²C positions in the presence of ¹²C citrate. Raman-SIP with D₂O labelling has been employed to probe metabolic activity of single cells without the need of cell replication. Our results show that *E. coli* alone in minimal medium with citrate as the sole carbon source had no metabolic activity, but became metabolically active in the presence of ADP1. Mass spectrometry-based metabolite footprint analysis suggests that putrescine and phenylalanine excreted by ADP1 cells may support the metabolic activity of *E. coli*. This study demonstrates that Raman-SIP with reverse labelling would be useful tool to probe metabolism of any carbon substrate, overcoming limitations when stable isotopic substrates are not readily available. Raman-SIP with D₂O labelling is a sensitive and reliable approach, which distinguishes metabolically active cells but not quiescent cells. This novel approach extends the application of Raman-SIP and demonstrates its potential application as a powerful tool for probing cellular metabolism and their metabolic interactions at the single cell level in microbial communities.

Table S1. List of the metabolites detected and identified by footprint GC-MS analysis

Name of compound	RI	RT	MSI ID level	ChEBI code
Tridecane	958.02	332.728	1	35998
Myristoleic acid	968.91	336.428	1	27781
Undecane	981.85	340.828	1	46342
Pyridine	1081.7	374.778	1	-
Decane	1201	415.328	1	41808
Oxalic acid	1231.8	425.828	2	16995
Glycerol	1284.5	443.728	1	17754
Phosphate	1418.9	489.428	1	18367
Unknown	1500.5	517.228	4	-
Unknown	1500.9	517.428	4	-
Unknown	1535.8	534.078	4	-
Cystathionine	1619.3	573.978	3	17755
Putrescine	1649	588.128	1	17148
Phenylalanine	1730.6	622.328	1	28044
Citric acid	1842.4	680.478	1	30769
Unknown	1901.9	708.528	4	-




Figure S1. Illustration of the reverse and multiple stable isotope labeling processes




Figure S2 (A). SCRS of *A. baylyi* ADP1 pure culture cells grown in MMC with/without 40% D₂O. After 3 h of incubation, there was no observable change of ¹²C-SCRS positions, and the C-D band appeared around 2187 cm⁻¹ in the presence of D₂O. **S2 (B).** SCRS of ¹³C-labeled *E. coli* DH5α-GFP pure culture cells grown in MMC with/without 40% D₂O. After 7 h of incubation, there was no observable change of ¹²C-SCRS positions and the C-D band did not appear around 2187 cm⁻¹ in the presence of D₂O.




Figure S3. Main features of Raman spectra of *A. baylyi* ADP1 and *E. coli* DH5 α -GFP pure cultures under no labelling (A), D₂O labelling (B), ¹³C labelling (C) and co-labelling of D₂O and ¹³C (D).




Figure S4. Growth curves of *A. baylyi* ADP1 and *E. coli* DH5α-GFP in MMC and MM medium under aerobic conditions. *A. baylyi* ADP1 could grow by using citrate as the sole carbon source but *E. coli* DH5α-GFP was unable to use citrate.




Figure S5. Combined image of the co-culture of DH5α-GFP & ADP1 cells




Figure S6. Cell counting results of *A. baylyi* ADP1 and *E. coli* DH5a-GFP in MMC with 40% D₂O. *A. baylyi* ADP1 could grow by using citrate as the sole carbon source. *E. coli* DH5a-GFP was able to be labeled by D₂O, but did not grow obviously.





Figure S7. (A) Growth curve of *E. coli* DH5a-GFP under different concentrations of D₂O. (B) Growth curve of *A. baylyi* ADP1 under different concentrations of D₂O.